

## Functional Electrical Stimulation System

### Background to the Invention

#### 5 1. Field of the Invention

The present invention relates to Functional Electrical Stimulation Systems.

#### 2. Background Information

10 Functional Electrical Stimulation (FES) systems artificially stimulate the muscles, and muscle groups, of persons through the use of electrical current in order to stimulate movement. As early as 1971 Liberson applied electrical stimulation to assist walking in patients with  
15 foot drop.

Most Functional Electrical Stimulation systems are designed for rehabilitation of spinal core injured or stroke patients. There are a number of FES systems on the market  
20 to assist foot drop in persons who have had a stroke. Typically, a heel switch is provided on the foot of a patient to indicate when the patient lifts their foot off the ground so that a controller can stimulate appropriate muscles to contract to raise the foot during the step.

Current systems suffer two main drawbacks. Firstly, it is not convenient to tune the electrical stimulation pulse characteristics to suit different patients during treatment because known apparatus do not have external controls to  
5 avoid unintentional tuning by patients. Secondly, known devices do not provide any information on a patient's adaptability to the device and improvement in step performance during the treatment.

#### 10 Summary of the Invention

It is an object of the present invention to provide a Functional Electrical Stimulation systems which assists with rehabilitation of stroke patients, or which at least  
15 offers a useful choice.

According to a first aspect of the invention there is provided an electrical stimulation device comprising: a sensor for detecting a movement event of a body part, an  
20 electrode for making electrical contact with an area of the body part, and a controller coupled to the sensor and electrode for receiving a sensor signal indicating the movement event, and for outputting to the electrode an output comprising a rise signal, a stimulation signal and a

fall signal, and programmed to record a duration of use and a number of movement events during the duration of use.

According to a second aspect of the invention there is  
5 provided an electrical stimulation device for controlling  
the movement of a body part comprising: a sensor for  
detecting a movement event of a body part, an electrode for  
making electrical contact with an area of the body part and  
for stimulating a muscle of the body part, a housing to be  
10 worn by a user of the device, receiver on the housing for  
receiving wireless signals from a remote unit, and a  
controller provided in the housing and coupled to the  
receiver for receiving stimulation data from the remote  
unit and storing the stimulation data in a stimulation  
15 file, and coupled to the sensor for receiving a sensor  
signal indicating the movement event, and for generating a  
control signal using the stimulation file in response to  
the movement event, and for outputting the control signal  
to the electrode.

20

Preferably, the controller is also programmed for  
generating a log file storing a duration of use and a  
number of movement events during the duration.

Preferably, the stimulation data includes a stimulation intensity level, a rise time, a stimulation time, and a fall time.

5 Preferably, the stimulation data also includes a pulse form, a triggering period, a triggering method and triggering criteria.

Preferably, the device further includes a computer  
10 removably coupled to the controller for downloading the stimulation file and log file, and for updating the stimulation file, and programmed to store data from the stimulation and log files in a database, and for outputting for display the stimulation data and the duration of use  
15 and the number of movement events.

Preferably, the computer is a Personal Digital Assistant.

Preferably, the database also includes information about  
20 the user of the device, and the computer is programmed for accessing the database by a Windows<sup>TM</sup> graphical user interface.

Preferably, the remote unit is a handheld remote control  
25 unit.

Preferably, the body part is a foot and the sensor is a heel switch.

5 Further aspects of the invention will become apparent from the following description, which is given by way of example only.

#### Brief Description of the Drawings

10

Embodiments of the invention will now be described with reference to the accompanying drawings in which:

Figure 1 is a schematic of the leg of a person with a  
15 Functional Electrical Stimulation (FES) device fitted,

Figure 2 is an exploded front and back view of a wearable control module for the FES device,

20 Figure 3 is an electrical schematic for the control module,

Figure 4 is a Graphical User Interface (GUI) for displaying and setting simulation data on a personal computer,

Figure 5 is a relationship diagram of database records for use in a Functional Electrical Stimulation rehabilitation system,

5 Figures 6 to 28 illustrate a Window<sup>TM</sup> GUI for assessing and updating the FES systems database,

Figure 29 is a perspective view of a Personal Digital Assistant (PDA) for interfacing with the FES system, and

10

Figure 30 is a handheld remote control unit for setting and updating simulation data of the FES device.

#### Description of the Preferred Embodiments

15

The preferred embodiment of the invention is an electrical stimulator device to assist in overcoming "foot drop" affecting a person who has had a stroke. A heel pressure switch is used to control the application of electrical  
20 pulses which, in turn, stimulate appropriate muscles to contract to raise the foot during step.

Referring to Figures 1 and 2, an electrical stimulation device for assisting with "foot drop" in a person includes  
25 a control module 2 which can be worn on the belt of a

patient. An input sensor in the form of heel switch 3 is provided to detect when the patient lifts their foot off the ground. Heel switch 3 is connected to control module 2 by wires 4. Electrodes 5 are positioned on the front part of the patient's leg and connected to control module 2 via wiring 4. The control module 2 includes a micro controller which is coupled to sensor 3 and electrode 5.

The control module 2 worn by a user includes a housing containing a micro controller for receiving an event signal from the heel switch 3 indicating that the user has lifted their foot. The controller generates a simulation signal and outputs it to the electrode 5 to simulate muscles which contract to lift the foot during the step. The housing includes a battery 6 for powering the controller. On a front portion of the housing are a battery light emitting diode 7 and a simulation light emitting diode 8 to indicate when a simulation is occurring. A test button 9 is provided for simulating the input from the heel switch 3 to cause a simulation output to the electrode 5. On one end of the housing is a rotary knob 10 for turning the module 2 on/off and adjusting the simulation intensity level. On the back of the housing is a belt clip adaptor 11 so that the housing can be worn by the user of the device.

Referring to Figure 3, a schematic diagram of the controller and associated electrical circuit is shown.

The computer also has a data logging function to record the  
5 hours of use and number of steps walked by the patient.  
This information is recorded in a log file which can be  
downloaded to a computer database or Personal Digital  
Assistant (PDA), shown in Figure 29. The information can  
be displayed on the computer or PDA screen.

10

By logging duration of use and number of steps taken the  
physician is able to know more about the daily life pattern  
of patients and so can help prescribe a more suitable time  
schedule of walking exercise for each individual user. The  
15 logged information can also be used to monitor overuse or  
insufficient walking exercise by patients. Patients can  
adapt to the device more effectively so that they can  
improve their quality of life by merging into the society  
again.

20

The housing also has a wireless communications receiver,  
such as an infra-red receiver, for receiving stimulation  
data from a remote handheld control unit, shown in Figure  
30. The remote control allows a physician to set  
25 stimulation parameters to fit individual patients. The



stimulation parameter data is sent to the controller and stored in a stimulation file which the controller uses to generate the stimulation signal when a step event is recorded by the heel switch. The wireless remote control  
5 makes stimulation parameter setting more effective and efficient. The physician can watch a patient and evaluate their walking characteristics while using the device. Immediate changes to the stimulation parameters can be made and the results observed without stopping the patient. As  
10 a consequence, patients will find it easier to adapt to the device and physicians will take less time in adjusting the device's parameters to fit different patients.

The computer and PDA software allows logged data to be  
15 retrieved from the device and recorded in a database for immediate or later analysis, and displayed on a screen. The PC software design is based on the familiar and user-friendly Windows™ Graphical User Interface. The Database contains a full range of patient data, as indicated in  
20 Figure 5, to allow the physician to retrieve and save all the patient records, which include all the parameters stored inside the stimulation devices.

Referring to Figures 6 to 28, Functional Electrical System  
25 (FES) Recording System comprises the computer database and

Windows GUI to allow storage of patient details and FES system parameters, and to allow monitoring patient's mobility as treatment progress.

5 Figure 6 illustrates the main window of the FES record system. The upper half of the window shows patient information including, for example, the patient's name, diagnosis, and address. The lower half of the window can be used to show other patient information arranged under  
10 seven pages accessible by tabs 125, 126, 127, 128, 129 130 and 131.

General Records 125 summaries important record information including appointment status, whether the patient is using  
15 the FES system or not and FES accessories bought so far. It provides a fast overview of the patient's information.

Appointment 126 shows information of appointments so that clinician can check progress of treatment.

20

FES Setting 127 shows the setting of the FES system during each patient visit to allow a continuous service to the patient by tracking the effect of different FES settings on the patient.

25

Product Lot 128 shows the date and number of accessories the patient has bought.

Photo 129 allows the clinician to upload digital images  
5 into the FES Recording System to supplement the patient's records visually. For example, an image may show the patient's leg with electrodes attached to show electrode attachment site.

10 Patient Status 130 shows the patient's physical information such as ROM, walking gait characteristics and muscle power.

Step Record 131 shows the walking pattern of the patient when using the FES System. The record includes the total  
15 number of steps, walking time (hours) and cadence per hour, per day or in the whole period.

Figure 7 shows a pop-up window used to add a new patient record. The pop-up window is revealed by clicking "Add  
20 New" 121 in the main window. The new patient information is added in each field 106 to 120. In order to saving time patient "Sex" 110 is preset to "Male" and "Hemiplegic Side" 117 is preset to "Right". To change the former to "Female" or the latter to "Left" or "Bilateral", click on arrow of  
25 the drop-down box.

For the box of "Injure Date" (118), you can input the date of injury in the form of "day/month/year" (d/m/y) or "month/year" (m/y), for those who forget the exact day. For  
5 patients who forget the date of injury, you can just choose delete (119) to not input this field. "Note" (120) is the field allow physician to record patient's history other than item provided.

10 You can edit a patient's information by click the icon "Edit" (122).

Patient records can be reviewed one by one by clicking "Previous" (123) or "Next" (124). To find a patient by  
15 name or patient ID, see "Search" screen below.

Referring to Figure 8, general records includes a patient latest appointment, use of the FES System or not, FESID number, the injured body segment, and accessories bought.

20

Referring to Figures 9 and 10, the FES Recording System allows easy and clear appointment booking. Click the "Add New" (134) and then select physician name (136, 137), date (138) and time of the appointment (139). Also the aim  
25 (140) for the appointment, whether it is a first visit or

follow up case, etc. Click "Edit" (135) to change an appointment.

Referring to Figures 11 and 12, you can review past  
5 settings of the FES system by highlighting the date of  
visit and click "View" (142). To set new settings for FES  
System or to change settings, select "Add New" (141) and  
input the new settings so that you can retrieve it anytime  
as shown in Figure 13.

10

Referring to Figure 14, physicians can record the number of  
accessories the patient has consumed by clicking "Add New"  
(143) or directly clicking on the specific icon to increase  
the item number. To decrease the number on a specific  
15 item, click the downward arrow under the corresponding icon  
directly and then click "Submit" (144).

Referring to Figure 15, physicians can upload photos to the  
FES Recording System. Photo can be, for example, electrode  
20 sites, skin condition before and after having electrode on  
patient's skin, toe clearance during walking with and  
without the application of FES on patient's affected foot.

Referring to Figure 16 and 17, physicians can update notes  
25 on a patient's physical condition during each visit. To

review previous patient's status, select the date of visit you would like to review and click "View" (146). "Patient Status Form" will be shown as in Figure 12. To update patient's physical status, click "Add New" (145), a blank  
5 Patient Status Form will then be shown. Select by directly clicking on the desired box. By clicking the icon (148), the physicians can type any other related information not included in the Form.

10 Referring to Figure 18, physicians can view a patient's walking pattern at a particular stage of treatment by highlighting the desired date of visit and click "View" (149).

15 Referring to Figure 19, 20 and 21, a patient's walking pattern using the FES System can be recorded. The number of steps and walking time in a selected hour, day or from a period of dates are shown in both numerical and graphic form. This information is logged by the wearable control  
20 module and downloaded to the database. Cadence is calculated by the software automatically. The physician can adjust the FES setting after monitoring patient's walking habit or after knowing patient's compliance.

Referring to Figure 22, physicians can add the name or location of their treatment facility by clicking "Add New" (151) and it will be showed on the front page of the system, as indicated by "Location" (104) in Figure 1. To  
5 change or delete, physicians can click "Edit" (152) or "Delete" (153) respectively.

Referring to Figure 23, physicians can add their name and their related information for booking appointment purpose  
10 by clicking "Add New" (155). To modify the information you can click "Edit" (156) or "Delete" (157) to delete it. FES Physician ID (159) will be created by the FES Recording System automatically.

15 Referring to Figure 24, you can add different diagnosis in the list by clicking "Add New" (164). To change or delete items in the list, you can click "Edit" (165) or "Delete" (166) respectively.

20 Referring to Figure 25, the balance of accessories in your center can be recorded. You can add or delete number of items by clicking "Add New" (168) or "Edit" (169) respectively.

Referring to Figure 26, you can search for a patient by using their PatientID (172), ID Card No. (173) or Name (174) and clicking "Search" (175).

5 Referring to Figure 27, "Filter" is used to select a specific group of patients with common criteria. For example the number of male patients aged 50 or more affected on his left side after a stroke.

10 Referring to Figure 28, "statistics" will tell you the status of the FES system. You can check the total number of FES systems sold or lent to patients in a certain period of time. You can also find out the number of FES system on loan to a patient. For example, to find out the number of  
15 FES Systems sold to patients from 28/02/2001 to 31/03/2003. Input "28/02/2001" in (177) and 31/03/2003 in (178) and then click "Search" (176). The result shown is 54 units (179).

20 Where in the foregoing description reference has been made to integers or elements having known equivalents then such are included as if individually set forth herein.

Embodiments of the invention have been described, however  
25 it is understood that variations, improvements or



modifications can take place without departure from the spirit of the invention or scope of the appended claims.